



## Data Sheet

# RETINOL-BINDING PROTEIN 4 (RBP4) (aa 169-183)

## ANTIBODY, POLYCLONAL

|                            |  |
|----------------------------|--|
| <b>Catalog no.:</b>        | A 1054.1 / A 1054.2  |
| <b>Immunogen:</b>          | Synthetic human RBP4 (aa 169-183) KLH-conjugated<br>(VHNGYCDGRSERNLL)  |
| <b>Synonyms:</b>           | Plasma retinol-binding protein (PRBP)  |
| <b>Swiss-Prot No:</b>      | P02753   |
| <b>Gene Information:</b>   | Gene Name: RBP4<br>Gene ID: 5950   |
| <b>Host:</b>               | Rabbit   |
| <b>Matrix:</b>             | Serum  |
| <b>Specificity:</b>        | Human RBP 4 (aa 169-183)   |
| <b>Contents:</b>           | 20 µl / 100 µl (lyophilized)<br>Resuspend in 20 µl / 100 µl aqua bidest.   |
| <b>Known applications:</b> | ELISA (1:2000)<br><br>This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls. |
| <b>Store at:</b>           | 2-8 °C (lyophilized); - 20 °C (dissolved)<br>Repeated thawing and freezing must be avoided   |

**Last updated on:** 16 March 2022

### For research use only

**Publishing research using A 1054? Please let us know so that we can cite your publication as a reference.**

