



Data Sheet

RETINOL-BINDING PROTEIN 4 (RBP4) (aa 169-179)

ANTIBODY, POLYCLONAL

Catalog no.: A 1057.1 / A 1057.2
Immunogen: Synthetic human RBP4 (aa 169-179) KLH-conjugated (VHNGYCDGRSE)
Synonyms: Plasma retinol-binding protein (PRBP)
Swiss-Prot No: P02753
Gene Information: Gene Name: RBP4

Gene ID: 5950

Host: Rabbit

Matrix: Serum

Specificity: Human RBP4 (aa 169-179)

Contents: 20 µl / 100 µl (lyophilized)
Resuspend in 20 µl / 100 µl aqua bidest.

Known applications: ELISA (1:2000)

This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.

Store at: 2-8 °C (lyophilized); - 20 °C (dissolved)
Repeated thawing and freezing must be avoided

Last updated on: 16 March 2022

For research use only

Publishing research using A 1057? Please let us know so that we can cite your publication as a reference.

