



## Data Sheet

# RETINOL-BINDING PROTEIN 4 (RBP4) (aa 169-181)

## ANTIBODY, POLYCLONAL

**Catalog no.:** A 1056.1 / A 1056.2  
**Immunogen:** Synthetic human RBP4 (aa 169-181) KLH-conjugated (VHNGYCDGRSERN)  
**Synonyms:** Plasma retinol-binding protein (PRBP)  
**Swiss-Prot No:** P02753  
**Gene Information:** Gene Name: RBP4

Gene ID: 5950

**Host:** Rabbit

**Matrix:** Serum

**Specificity:** Human RBP4 (aa 169-181)

**Contents:** 20 µl / 100 µl (lyophilized)  
Resuspend in 20 µl / 100 µl aqua bidest.

**Known applications:** ELISA (1:2000)

This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.

**Store at:** 2-8 °C (lyophilized); -20 °C (dissolved)  
Repeated thawing and freezing must be avoided

**Last updated on:** 16 March 2022

### For research use only

**Publishing research using A 1056? Please let us know so that we can cite your publication as a reference.**

