



## Data Sheet

# RETINOL-BINDING PROTEIN 4 (RBP4) (aa 169-182)

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	A 1055.1 / A 1055.2
<b>Immunogen:</b>	Synthetic human RBP4 (aa 169-182) KLH-conjugated (VHNGYCDGRSERNL)
<b>Synonyms:</b>	Plasma retinol-binding protein (PRBP)
<b>Swiss-Prot No:</b>	P02753
<b>Gene Information:</b>	Gene Name: RBP4 Gene ID: 5950
<b>Host:</b>	Rabbit
<b>Matrix:</b>	Serum
<b>Specificity:</b>	Human RBP4 (aa 169-182)
<b>Contents:</b>	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
<b>Known applications:</b>	ELISA (1:2000) <small>This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.</small>
<b>Store at:</b>	2-8 °C (lyophilized); -20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	16 March 2022

### For research use only

**Publishing research using A 1055? Please let us know so that we can cite your publication as a reference.**

