



## Data Sheet

# HUMAN RETINOL-BINDING PROTEIN 4 (RBP4)

## NATIVE

<b>Catalog no.:</b>	A 1059AG.1 / A 1059AG.2
<b>Name:</b>	Human retinol-binding protein 4 (RBP4)
<b>Synonyms:</b>	Plasma retinol-binding protein (PRBP)
<b>MW:</b>	21 kDa
<b>Swiss-Prot No:</b>	P02753
<b>Gene Information:</b>	Gene Name: RBP4 GeneID: 5950
<b>Source:</b>	Urine from patients with renal tubular proteinuria
<b>Storage buffer:</b>	20 mM NH <sub>4</sub> CO <sub>3</sub>
<b>Purity:</b>	> 98%
<b>Contents:</b>	100 µg / 1 mg (lyophilized) Resuspend in phosphate buffer, pH >7.0 containing 15 mM NaCl for required concentration
<b>Store at:</b>	- 20 °C Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	04 May 2022

### For research use only

**Publishing research using A 1059AG? Please let us know so that we can cite your publication as a reference.**

