

## Data Sheet

# RISANKIZUMAB

## ANTIBODY, RECOMBINANT MONOCLONAL, COUPLED TO HSA

<b>Catalog no.:</b>	A 1235.3 / A 1235.4
<b>Immunogen:</b>	Risankizumab
<b>Synonyms:</b>	Skyrizi
<b>DrugBank:</b>	DB14762
<b>Host:</b>	Human cell line
<b>Clone no.:</b>	YU1220-E02
<b>Form:</b>	Monovalent single chain variable fragment (scFv) coupled to HSA
<b>Matrix:</b>	PBS pH 7.4
<b>Specificity:</b>	Risankizumab variable domain
<b>Contents:</b>	10 µg / 100 µg (liquid)
<b>Known applications:</b>	ELISA (capture) <small>This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.</small>
<b>Store at:</b>	-20 °C
<b>Last updated on:</b>	31 October 2024

### For research use only

**Publishing research using A 1235? Please let us know so that we can cite your publication as a reference.**

