



Data Sheet

HUMAN SECRETIN (aa 1-13)

ANTIBODY, POLYCLONAL

Catalog no.:	AE1038.1 / AE1038.2
Immunogen:	Synthetic human Secretin (aa 1-13) poly-Lysin conjugated (HSDGTFTSELSRL)
Swiss-Prot No:	P09683
Gene Information:	Gene Name: SCT GenelD: 6343
Host:	Rabbit
Matrix:	Serum
Specificity:	Synthetic human Secretin (aa 1-13)
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA (1:5000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	20 April 2022

For research use only

Publishing research using AE1038? Please let us know so that we can cite your publication as a reference.

