

## Data Sheet

# HUMAN SUBSTANCE P

## SYNTHETIC

<b>Catalog no.:</b>	A 6103AG.1 / A 6103AG.2
<b>Name:</b>	Human Substance P
<b>Synonyms:</b>	Protachykinin-1 (aa 58-68)
<b>Sequence:</b>	RPKPQQFFGLM
<b>MW:</b>	1.4 kDa
<b>Swiss-Prot No:</b>	P20366
<b>Gene Information:</b>	Gene Name: TAC1, NKA, NKNA, or TAC2 GeneID: 6863
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	100 µg/ 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
<b>Store at:</b>	- 20 °C Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	06 May 2022

**For research use only**

**Publishing research using A 6103AG? Please let us know so that we can cite your publication as a reference.**