

Data Sheet

HUMAN THYMOPENTIN SYNTHETIC

Catalog no.:	A 9590AG.1 / A 95901AG.2
Name:	Human Thymopentin (derivative of thymopoietin)
Synonyms:	TP5
Sequence:	RKDVY
MW:	0.68 kDa
Swiss-Prot No:	P42166
Gene Information:	Gene Name: TMPO Gene ID: 7112
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration.
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	11 May 2022

For research use only

Publishing research using A 9590AG? Please let us know so that we can cite your publication as a reference.

