

## Data Sheet

# HUMAN PROTHYMOSIN $\alpha$ (aa 103-111) SYNTHETIC

<b>Catalog no.:</b>	A 9571AG.1 / A 9571AG.2
<b>Name:</b>	Human Prothymosin $\alpha$ (aa 103-111)
<b>Sequence:</b>	KKQKTDEDD
<b>MW:</b>	1.11 kDa
<b>Swiss-Prot No:</b>	P06454
<b>Gene Information:</b>	Gene Name: PTMA, TMSA Gene ID: 5757
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	100 $\mu$ g / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration.
<b>Store at:</b>	- 20 °C Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	4 January 2023

**For research use only**

**Publishing research using A 9571AG? Please let us know so that we can cite your publication as a reference.**

