

Data Sheet

HUMAN PROTHYMOSIN α_1

SYNTHETIC

Catalog no.:	A 9541AG.1 / A 9541AG.2
Name:	Human ProThymosin α_1 (aa 2-111) (AcSDAAVDTSSSEITTKDLKEKKEVVEEAENGRDAPANGNAENEENGQEADNEV DEEEEEGGEEEEEEEEEGDGEEDGDEDEEAESATGKRAAEDDEDDVDTKKQKT DEDD)
MW:	12.05 kDa
Swiss-Prot No:	P06454
Gene Information:	Gene Name: PTMA GenelD: 5757
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	10 μ g/100 μ g (lyophilized) Resuspend in aqua bidest for required concentration
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	11 May 2022

For research use only

Publishing research using A 9541AG? Please let us know so that we can cite your publication as a reference.