

Data Sheet

HUMAN THYMOSIN α 1

SYNTHETIC

Catalog No.:	A 9510AG.1 / A 9510AG.2
Name:	Human Thymosin α 1 (aa 2-29)
Synonyms:	Prothymosin α (aa 2-29)
Sequence:	AcSDAAVDTSSSEITTKDLKEKKEVVEEAEN
MW:	3.1 kDa
Swiss-Prot No:	P06454
Gene Information:	Gene Name: PTMA, TMSA Gene ID: 5757
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 μ g / 1 mg (lyophilized) Resuspend in in 0.1% acetic acid for required concentration
Store at:	-20°C Repeated thawing and freezing must be avoided
Last updated on:	11 May 2022

For research use only

Publishing research using A 9510AG? Please let us know so that we can cite your publication as a reference.

