

Data Sheet

HUMAN THYMOSIN β 10 (aa 13-42) SYNTHETIC

| | |
|--------------------------|---|
| Catalog no.: | A 9568AG.1 / A 9568AG.2 |
| Name: | Human Thymosin β 10 (aa 13-42) |
| Sequence: | DKAKLKKTTETQEKNTLPTKETIEQEKRSEI |
| MW: | 3.55 kDa |
| Swiss-Prot No: | P63313 |
| Gene Information: | Gene Name: TMSB10, PTMB10, THYB10 Gene ID: 9168 |
| Source: | Synthetic |
| Purity: | > 95% (HPLC) |
| Contents: | 100 μ g / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration. |
| Store at: | - 20 °C Repeated thawing and freezing must be avoided |
| Last updated on: | 11 May 2022 |

For research use only

Publishing research using A 9568AG? Please let us know so that we can cite your publication as a reference.

