

Data Sheet

HUMAN THYMOSIN β 10 SYNTHETIC

Catalog no.:	A 9563AG.1 / A 9563AG.2
Name:	Human Thymosin β 10
Sequence:	AcADKPDMGGEIASFDKAKLKKTTETQEKNLPTKETIEQEKRSSEIS
MW:	4.9 kDa
Swiss-Prot No:	P63313
Gene Information:	Gene Name: TMSB10, PTMB10, THYB10 Gene ID: 9168
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 μ g / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration.
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	11 May 2022

For research use only

Publishing research using A 9563AG? Please let us know so that we can cite your publication as a reference.

