



## Data Sheet

# HUMAN THYMOSIN $\beta$ 10 (aa 13-43) SYNTHETIC

<b>Catalog no.:</b>	A 9565AG.1 / A 9565AG.2
<b>Name:</b>	Human Thymosin $\beta$ 10 (aa 13-43)
<b>Sequence:</b>	DKAKLKKTTETQEKNTLPTKETIEQEKRSSEIS
<b>MW:</b>	3.63 kDa
<b>Swiss-Prot No:</b>	P63313
<b>Gene Information:</b>	Gene Name: TMSB10, PTMB10, THYB10 Gene ID: 9168
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	100 $\mu$ g / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration.
<b>Store at:</b>	- 20 °C Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	11 May 2022

**For research use only**

**Publishing research using A 9565AG? Please let us know so that we can cite your publication as a reference.**

