

Data Sheet

HUMAN THYMOSIN β 15 (aa 1-45)

SYNTHETIC

Catalog no.:	A 9575AG.1 / A 9575AG.2
Name:	Human Thymosin β 15 (aa 1-45)
Synonyms:	NB thymosin beta, Thymosin-like protein 8
Sequence:	MGDRPDLSEVERFDKSKLKKTTITEVKNTLPSKETIEQEKEFVKRS
MW:	5.3 kDa
Swiss-Prot No:	P0CG34
Gene Information:	Gene Name: TMSB15A, TMSL8, TMSNB Gene ID: 11013, 286527
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 μ g / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration.
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	11 May 2022

For research use only

Publishing research using A 9575AG? Please let us know so that we can cite your publication as a reference.

