



Data Sheet

HUMAN THYMOSIN β 4 (aa 38-43)

SYNTHETIC

Catalog no.:	A 9550AG.1 / A 9550AG.2
Name:	Human Thymosin β 4 (aa 38-43)
Sequence:	KQAGES
Synonyms:	T beta-4, Fx
MW:	0.62 kDa
Swiss-Prot No:	P62328
Gene Information:	Gene Name: TMSB4X, TB4X, THYB4, TMSB4 Gene ID: 7114
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 μ g / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	11 May 2022

For research use only

Publishing research using A 9550AG? Please let us know so that we can cite your publication as a reference.

