

Data Sheet

HUMAN THYMOSIN β 4 (C-TERMINAL BIOTINYLATED)

SYNTHETIC

Catalog No.:	A 9525AG.1 / A 9525AG.2
Name:	Human Thymosin β 4 (C-terminal biotinylated)
Synonyms:	T beta-4, Fx
Sequence:	SDKPDMAEIEKFDKSKLKKTTETQEKNPLPSKETIEQEKGAGES- <u>L</u> -Biotin Underlined amino acids were added or modified for technical reasons.
MW:	5.2 kDa
Swiss-Prot No:	P62328
Gene Information:	Gene Name: TMSB4X, TB4X, THYB4, TMSB4 GeneID: 7114
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 μ g / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
Store at:	-20°C Repeated thawing and freezing must be avoided
Last updated on:	11 May 2022

For research use only

Publishing research using A 9525AG? Please let us know so that we can cite your publication as a reference.

