

## Data Sheet

# HUMAN THYMOSIN $\beta$ 4 (aa 1-11)

## SYNTHETIC

<b>Catalog no.:</b>	A 9523AG.1 / A 9523AG.2
<b>Name:</b>	Human Thymosin $\beta$ 4 (aa 1-11)
<b>Sequence:</b>	AcSDKPDMAEIEK
<b>Synonyms:</b>	T beta-4, Fx
<b>MW:</b>	1.26 kDa
<b>Swiss-Prot No:</b>	P62328
<b>Gene Information:</b>	Gene Name: TMSB4X, TB4X, THYB4, TMSB4 Gene ID: 7114
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	100 $\mu$ g / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
<b>Store at:</b>	- 20 $^{\circ}$ C Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	11 May 2022

**For research use only**

**Publishing research using A 9523AG? Please let us know so that we can cite your publication as a reference.**

