

Data Sheet

HUMAN THYROGLOBULIN (TG) NATIVE

Catalog no.:	AK3122AG.1 / AK3122AG.2
Name:	Human native Thyroglobulin
MW:	660 kDa (dimer)
Swiss-Prot No:	P01266
Gene Information:	Gene Name: TG Gene ID: 7038
Source:	Human tissue
Storage buffer:	20 mM NH ₄ CO ₃
Purity:	> 95% (SDS-PAGE)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in PBS, 150 mM NaCl pH 7.0 for required concentration.
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	20 May 2022

For research use only

Publishing research using AK3122AG? Please let us know so that we can cite your publication as a reference.

