



Data Sheet

HUMAN THYROXINE (T4)

ANTIBODY, POLYCLONAL

| | |
|----------------------------|--|
| Catalog no.: | A 1037.1 / A 1037.2 |
| Immunogen: | Synthetic T4 BSA-conjugated |
| Host: | Rabbit |
| Matrix: | Serum |
| Specificity: | T4 |
| Contents: | 20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest. |
| Known applications: | RIA (1:1000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls. |
| Store at: | 2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided |
| Last updated on: | 16 March 2022 |

For research use only

Publishing research using A 1037? Please let us know so that we can cite your publication as a reference.

