



Data Sheet

HUMAN TIP 39

ANTIBODY, POLYCLONAL

Catalog no.:	A 1051.1 / A 1051.2
Immunogen:	Synthetic human TIP 39 (aa 62-100) bTG conjugated (SLALADDAAFRERARLLAALERRHWLNSYMHKLLVLDAP)
Synonyms:	Tuberoinfundibular peptide of 39 residues, Parathyroid hormone 2
Swiss-Prot No:	Q96A98
Gene Information:	Gene Name: PTH2, TIP39, TIPF39 Gene ID: 113091
Host:	Rabbit
Matrix:	Serum
Specificity:	Human, bovine, and canine TIP 39
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA (1:4000); RIA (1: 2000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	16 March 2022

For research use only

Publishing research using A 1051? Please let us know so that we can cite your publication as a reference.

