



Data Sheet

HUMAN TIP 39

SYNTHETIC

Catalog no.:	A 1050AG.1 / A 1050AG.2
Name:	Synthetic human TIP 39 (aa 62-100)
Sequence:	SLALADDAAFRERARLLAALERRHWLNSYMHKLLVLDAP
Synonyms:	Tuberoinfundibular peptide of 39 residues, Parathyroid hormone 2
Swiss-Prot No:	Q96A98
Gene Information:	Gene Name: PTH2, TIP39, TIPF39 GeneID: 113091
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in aqua bidest. for required concentration
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Reference:	1. Hansen IA, Jakob O, Wortmann S, Arzberger T, Allolio B, Blind E (2002). Characterization of the human and mouse genes encoding the tuberoinfundibular peptide of 39 residues, a ligand of the parathyroid hormone receptor family. <i>J Endocrinol</i> 174 (1): 95-102.
Last updated on:	04 May 2022

For research use only

Publishing research using A 1050AG? Please let us know so that we can cite your publication as a reference.

