



Data Sheet

HUMAN TISSUE PLASMINOGEN ACTIVATOR (tPA)

ANTIBODY, POLYCLONAL

Catalog no.:	A 7641.1 / A 7641.2
Immunogen:	Purified human tissue plasminogen activator (tPA)
Swiss-Prot No:	P00750
Gene Information:	Gene Name: PLAT Gene ID: 5327
Host:	Goat
Matrix:	Serum
Specificity:	Free human tPA and tPA-inhibitor-complexes. There was no cross reactivity obtained with human uPA and PAI-2.
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA, immunohistochemistry (cryosections, 1:50-100; paraffin sections, 1:50) Recommended procedure for paraffin sections: Pretreatment with 0.05% Pronase E (Sigma) for 20 minutes. This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	07 April 2022

For research use only

Publishing research using A 7641? Please let us know so that we can cite your publication as a reference.

