



Data Sheet

HUMAN TRANSFORMING GROWTH FACTOR α [TGF α]

ANTIBODY, POLYCLONAL

Catalog no.: AS1022.1 / AS1022.2
Immunogen: Recombinant human Transforming Growth Factor α
Synonyms: TGF-alpha, EGF-like TGF (ETGF), TGF type 1
Swiss-Prot No: P01135

Gene Information: Gene Name: TGFA
GenelD: 7039

Host: Rabbit

Matrix: Serum

Specificity: Human Transforming Growth Factor α

Contents: 20 μ l / 100 μ l (lyophilized)
Resuspend in 20 μ l / 100 μ l aqua bidest.

Known applications: Western Blot

This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.

Store at: 2-8 °C (lyophilized); -20 °C (dissolved)
Repeated thawing and freezing must be avoided

Last updated on: 28 April 2022

For research use only

Publishing research using AS1022? Please let us know so that we can cite your publication as a reference.

