

Data Sheet

HUMAN TRANSFORMING GROWTH FACTOR β 1 (TGF β 1) RECOMBINANT

Catalog no.:	A 9600AG.3 / A 9600AG.4
Name:	Human TGF β 1 (aa 250-361) (ALDTNYCFSSTEKNCCVRQLYIDFRKDLGWKWIHEPKGYHANFCLGPCPYIWSL DTQYSKVLALYNQHNP GASAAPCCVPQALEPLPIVYVYVGRKPKVEQLSNMIVRSC KCS)
Swiss-Prot No:	P01137
Gene Information:	Gene Name: TGF β 1, TGF β GeneID: 7040
Source:	HEK293 cells
Storage buffer:	0.1% trifluoroacetic acid
Purity:	> 98% (SDS-PAGE and HPLC)
Contents:	50 μ g / 100 μ g (lyophilized) Resuspend in 10 mM citric acid pH 3.0 to a concentration of 0.1-1.0 mg/ml. Do not vortex. For extended storage dilute further in a buffer containing a carrier protein (for example 0.1% BSA).
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	18 September 2023

For research use only

Publishing research using A 9600AG? Please let us know so that we can cite your publication as a reference.

