

## Data Sheet

# HUMAN TISSUE TRANSGLUTAMINASE [TG2]

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	AK3010.1 / AK3010.2
<b>Immunogen:</b>	Recombinant human TissueTransglutaminase [TG2]
<b>Synonyms:</b>	Protein-glutamine gamma-glutamyltransferase 2, EC=2.3.2.13, Transglutaminase C, Transglutaminase H, Transglutaminase-2
<b>Swiss-Prot No:</b>	P21980
<b>Gene Information:</b>	Gene Name: TGM2 Gene ID: 7052
<b>Host:</b>	Rabbit
<b>Matrix:</b>	Serum
<b>Specificity:</b>	Human TissueTransglutaminase [TG2]
<b>Contents:</b>	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
<b>Known applications:</b>	ELISA (> 1:4000)  This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	06 October 2023

### For research use only

**Publishing research using AK3010? Please let us know so that we can cite your publication as a reference.**

