

Data Sheet

BOVINE TROPONIN T1 (HRP LABELED) (aa 16-31) SYNTHETIC

Catalog no.:	A 1002AG.1 / A 1002AG.2
Name:	Bovine Troponin T1 (aa 16-31)
Synonyms:	Troponin T, fast skeletal muscle (TnTf)
Sequence:	APPPAEVPEVHEEVH
MW:	1.64 kDa
Swiss-Prot No:	Q8MKI3
Gene Information:	Gene Name: Tnnt3 Gene ID: 282096
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration.
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	4 May 2022

For research use only

Publishing research using A 1002AG? Please let us know so that we can cite your publication as a reference.

