



Data Sheet

ANTI HUMAN UROTENSIN II

ANTIBODY, POLYCLONAL

Catalog no.:	A 8020.1 / A 8020.2
Immunogen:	Synthetic human Urotensin II, KLH-conjugated (ETPD CFWKYCV)
Synonyms:	Urotensin-2, U-II, UII
Swiss-Prot No:	O95399
Gene Information:	Gene Name: UTS2 GenelD: 10911
Host:	Rabbit
Matrix:	Serum
Specificity:	Synthetic human Urotensin II There was no cross reactivity obtained with human Endothelin (aa 1-21) (< 0.1%); human big Endothelin (aa 1-38) (< 0.1%); human Urodilatin; human ANP (1-28) (< 0.1%); human BNP (aa 1-32) (<0.1%); Adrenomedulin (aa 1-52) (< 0.1%); Relaxin (< 0.1%); or Somatostatin (< 0.1%).
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA (1:3000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	08 April 2022

For research use only

Publishing research using A 8020? Please let us know so that we can cite your publication as a reference.

