

## Data Sheet

# HUMAN UROTENSIN II

## SYNTHETIC

<b>Catalog no.:</b>	A 8024AG.1 / A 8024AG.2
<b>Name:</b>	Human Urotensin II, bTG conjugated
<b>Sequence:</b>	ETPDCFWKYCV
<b>Synonyms:</b>	Urotensin-2
<b>Swiss-Prot No:</b>	O95399
<b>Gene Information:</b>	Gene Name: UTS2 GeneID: 10911
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	10 µg / 100 µg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
<b>Store at:</b>	- 20 °C Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	06 May 2022

**For research use only**

**Publishing research using A 8024AG? Please let us know so that we can cite your publication as a reference.**

