



Data Sheet

HUMAN VASOACTIVE INTESTINAL PEPTIDE-RECEPTOR 2 [VIP-R2] (aa 3-14)

ANTIBODY, POLYCLONAL

Catalog no.:	AE6302.1 / AE6302.2
Immunogen:	Synthetic human VIP-R2 (aa 3-14) (RFHLEIQEEETK) polyLysine coupled
Synonyms:	Helodermin-preferring VIP receptor, Pituitary adenylate cyclase-activating polypeptide type III receptor (PACAP-R3), VPAC2
Swiss-Prot No:	P41587
Gene Information:	Gene Name: VIPR2, VIP2R GenelD: 7434
Host:	Rabbit
Matrix:	Serum
Specificity:	Human VIP-R2 (aa 3-14)
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Store at:	2-8 °C (lyophilized); -20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	22 April 2022

For research use only

Publishing research using AE6302? Please let us know so that we can cite your publication as a reference.

