



## Data Sheet

# VIBRIO CHOLERAE ZONA OCCLUDENS TOXIN (ZOT)

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	A 1225.1 / A 1225.2
<b>Immunogen:</b>	Synthetic vibrio cholerae Zot (aa 291-298) (GRLCVQDG) KLH-conjugated
<b>Synonyms:</b>	Zonular occludens toxin
<b>Swiss-Prot No:</b>	P38442
<b>Gene Information:</b>	Gene Name: zot GeneID: 2613964
<b>Host:</b>	Rabbit
<b>Matrix:</b>	Serum
<b>Specificity:</b>	Vibrio cholerae Zot (aa 291-298)
<b>Contents:</b>	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
<b>Known applications:</b>	ELISA (1:1000)  This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	16 March 2022

### For research use only

**Publishing research using A 1225? Please let us know so that we can cite your publication as a reference.**

