


Data Sheet

3-EPI-VITAMIN D₃

SYNTHETIC

Catalog no.: AD1006AG.2
Name: 3-Epi-Vitamin D₃
Synonyms: CAS 57651-82-8



Structure:

Molecular formula: C₂₇H₄₄O

MW: 384.65

Source: Synthetic

Purity: >98%

Appearance: Colourless, amorphous solid


Solubility: Soluble in EtOH/MeOH, DMF and DMSO, well soluble in halogenated hydrocarbons, not soluble in water

Contents: 1 mg (lyophilized)
Resuspend in 99% EtOH p.a. for required concentration

Store at: - 20 °C, protect from heat and light

Last updated on: 12 May 2022





For research use only

Publishing research using AD1006AG? Please let us know so that we can cite your publication as a reference.



Immundiagnostik AG

Stubenwald-Allee 8a · 64625 Bensheim · Germany

Phone: +49 6251 70190-0 · Fax: +49 6251 70190-363 · dept.immuochemicals@immundiagnostik.com · www.immundiagnostik.com