



Data Sheet

1 α -HYDROXYVITAMIN D₂

SYNTHETIC

Catalog no.: AD1016AG.2
Name: 1 α -Hydroxyvitamin D2
Synonyms: CAS 54573-75-0



Structure: 

Molecular formula: C₂₈H₄₄O₂

MW: 412.66

Source: Synthetic

Purity: >98%

Appearance: Colourless, amorphous solid

Solubility: Soluble in EtOH/MeOH, DMF and DMSO, well soluble in halogenated hydrocarbons, not soluble in water

Contents: 1 mg (lyophilized)
Resuspend in 99% EtOH p.a. for required concentration

Store at: - 20 °C, protect from heat and light

Last updated on: 24 May 2022

For research use only

Publishing research using AD1016AG? Please let us know so that we can cite your publication as a reference.

