


Data Sheet

3-EPI-25-HYDROXYVITAMIN-D₃

SYNTHETIC

Catalog no.:	AD1007AG.2
Name:	3-Epi-25-Hydroxyvitamin-D ₃
Synonyms:	CAS 73809-05-9



Structure:	
Molecular formula:	C ₂₇ H ₄₄ O ₂
MW:	400.65
Source:	Synthetic
Purity:	>98%
Appearance:	Colourless, amorphous solid
Solubility:	Soluble in EtOH/MeOH, DMF and DMSO, well soluble in halogenated hydrocarbons, not soluble in water
Contents:	1 mg (lyophilized) Resuspend in 99% EtOH p.a. for required concentration
Store at:	- 20 °C, protect from heat and light
Last updated on:	12 May 2022

For research use only

Publishing research using AD1007AG? Please let us know so that we can cite your publication as a reference.

