



Data Sheet

HUMAN VITAMIN D BINDING PROTEIN

ANTIBODY, POLYCLONAL, PEROXIDASE-LABELED

Catalog no.:	A 1020.5 / A 1020.6
Immunogen:	Human GC Globulin
Synonyms:	Vitamin D-binding protein (VDBP, DBP), Gc protein-derived macrophage activating factor1 (Gc-MAF, GcMAF1); Gc-globulin, Group-specific component (Gc1), Vitamin D-binding protein-macrophage activating factor (DBP-maf)
Swiss-Prot No:	P02774
Gene Information:	Gene Name: GC GenelD: 2638
Host:	Rabbit
Matrix:	IgG; purified from serum; peroxidase labeled
Specificity:	GC Globulin
Contents:	100 µg / 1 mg (liquid)
Known applications:	ELISA This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	- 20 °C
Last updated on:	04 April 2025

For research use only

Publishing research using A 1020? Please let us know so that we can cite your publication as a reference.

