

## Data Sheet

# HUMAN VITAMIN D BINDING PROTEIN NATIVE

<b>Catalog no.:</b>	A 1020AG.1 / A 1020AG.2
<b>Name:</b>	Vitamin D-binding protein
<b>Synonyms:</b>	Human GC Globulin (VDBP, DBP), Gc protein-derived macrophage activating factor1 (Gc-MAF, GcMAF1); Gc-globulin, Group-specific component (Gc1), Vitamin D-binding protein-macrophage activating factor (DBP-maf)
<b>MW:</b>	51.25 kDa
<b>Swiss-Prot No:</b>	P02774
<b>Gene Information:</b>	Gene Name: GC GenelD: 2638
<b>Source:</b>	Purified from human plasma
<b>Storage buffer:</b>	150 mM NaCl, 20 mM PBS pH 7.4
<b>Purity:</b>	≥95% by SDS-PAGE
<b>Contents:</b>	100 µg / 1 mg (lyophilized) Resuspend in aqua bidest. for required concentration.
<b>Store at:</b>	-20°C Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	04 May 2022

**For research use only**

**Publishing research using A 1020AG? Please let us know so that we can cite your publication as a reference.**

